

Sedens 400

Hybrid Seat Cushion

A light weight and compact seat cushion

The Sedens 400 aligns with the 2019 clinical guideline, offering cushions and back supports that prevent tissue injury and manage pressure effectively. Featuring a hybrid design with alternating modes, it redistributes pressure effectively. The innovative air valve allows control over immersion levels for optimal pressure distribution. Designed with careful selection of materials and covers, the Sedens 400 ensures consistent performance, reducing the risk of pressure injuries. Regular inspection and evaluation of our seating solutions guarantee long-term support and effectiveness.



Easy Setup Cushion

The cushions are designed to be easily assembled, providing instant comfort without the need for electricity.



Anti-Bottoming Design

The air cell-like structure foam encasement is designed to reduce the risk of bottoming out for patient safety.

Sedens 400

A cost-effective seat cushion solution that meets patients' clinical requirements



Continual Pressure Relief Care

The cushion can connect to pump and offers an alternating mode that effectively redistributes pressure.

Hybrid Cells

Single-piece pad design with six air cell-like structure and foam encasement.



Cap with Pressure Adjustment Valve

Innovative air valve to control the amount of immersion for optimal pressure distribution.

4 Way Stretch Top Cover

Offering low shear, low friction, vapor permeability, and moisture protection, coupled with optimal immersion and envelopment.



Non-Slip Base

Prevents from sliding out of position, giving peace of mind to the user.



Water Resistance



Flame Resistance



Biocompatibility



Vapour Permeable



Fungistatic



4-Way Stretch

Specifications

Model	Sedens 400				
Cushion	Dimension	Cells	Length	Width	Height
		6	17.3 in	16.1 in	3.5 in
	Weight	2.4 lb			
	Top cover material	Poly/PU			
	Cell material	Nylon/PU			
	Base material	PVC			
	Maximum patient weight	264 lb			
	HCPCS code	E2624			
	Ordering Information (For North America Only)				
Cushion	PM02101				

Cushion flame retardant standards (EN597-1, EN597-2), RoHS, WEEE
Specifications are subject to change for improvement without notice

Wellell America Corp.

927 Mariner Street,
Brea, CA 92821, USA
T +1 714 671 3818
F +1 714 494 8173
info.usa@wellell.com
www.wellell.com
Print-2024/08 All rights reserved