



# Sedens 400

## Hybrid Seat Cushion

- Alternating cushion works with Apex dynamic mattress pumps to provide effective seating solution
- Lightweight and portable, perfectly to be used in chairs
- Hybrid design with alternating mode to redistribute pressure effectively
- Innovative air valve to control the amount of immersion for optimal pressure distribution



Care for a Healthy Life

**APEX**

# Sedens 400

## Hybrid Seat Cushion



**24 Hours Pressure Area Care**  
Alternating cushion works with Apex dynamic mattress pumps to provide effective seating solution




**Hybrid Cells**  
6 air-cell-like one-piece pad with foam encased in to reduce risk of bottoming out



**Positioning Cushion**  
Innovative air valve to control the amount of immersion for optimal pressure distribution

**Non-Slip Base**  
Prevents from sliding out of position, giving peace of mind to the user

**4-Way Stretch Cover**  
Fluid resistant, low shear and vapor permeable top cover protect skin from friction and moisture



### Stretch Cover



Water Resistance



Flame Resistance



Biocompatibility



Vapour Permeable



Fungistatic



4-Way Stretch

### Specification

Cushion Dimension	17.3" x 16.1" x 3.5" / 44 x 41 x 9 cm
Cell Height	3.5"
Weight	2.4 lb / 1.11 Kg
Cover Material	4-Way Stretch
Cell Material	Nylon / PU
Base Material	Polyester PVC
Arrangement	6 Hybrid Cells (Cellular bladders filled with foam)
Flame Retardant Standards	EN597-1 & EN597-2

*Specifications are subject to change for improvement without notice*

### Ordering Information (For North America Only)

Model	Sedens 400
Cushion	PM02101

Distributed by

Apex Medical USA  
927 Mariner Street,  
Brea, CA 92821, USA  
T +1 714 671 3818  
F +1 714 494 8173  
sales.usa@apexmedicalcorp.com  
www.apexmedicalcorp.com  
Print-2019/05-500 All rights reserved