

- Lithium-ion battery of long duration that allows an autonomy up to 12 hours
- Visual and audible alarm enhances patients safety and security









# Sedens 500

## Seat Cushion with Battery Power



#### Specification

эрооннос	20011							
Pump	Dimension	9.8" x 4.7" x 2.2" / 25 x 12.2 x 5.7 cm						
	Weight	2.2 lb / 1.0 Kg						
	Case Material	Flame retardant ABS						
	Supply Voltage	12V DC, 1150 mA						
	Operating Cycle	10/15/20 min. selectable						
	Li-ion Battery Pack	3.78V Li-ion Battery x 2						
	Battery Operation Time	Approx. 12 hours						
	Battery Recharge Time	Approx. 4.5 hours						
Cushion	Dimension	16" x 16" x 4"	17" x 17" x 4"	18" x 16" x 4"				
		41 x 42 x 10.2cm	43 x 43 x 10.2cm	46 x 41 x 10.2cm				
		18" x 18" x 4"	20" x 16" x 4"	20" x 18" x 4"				
		46 x 46 x 10.2cm	51 x 41 x 10.2cm	51 x 46 x 10.2cm				
	Weight	2.4 lb / 1.1 Kg						
	Material	Cover: 4-Way Stretch; Cell: Nylon/PU & PU Foam; Base: Nylon/P						
	Maximum Patient Weight	260 lb / 120 kg						
	Flame Retardant Standards	EN597-1 & EN597-2						
HCPCS code		E2610						

 $Specifications \ are \ subject \ to \ change \ for \ improvement \ without \ notice$ 

#### Ordering Information (For North America Only)

Model	Sedens 500							
	16" x 16"	17" x 17"	18" x 16"	18" x 18"	20" x 16"	20" x 18"		
System Standard Cover	PD01905	PD01901	PD01902	PD01906	PD01903	PD01904		
Cushion Standard Cover	PM02901	PM02102	PM02701	PM02902	PM02801	PM02802		
Pump	PD01001							

Distributed by

### Apex Medical USA

927 Mariner Street, Brea, CA 92821, USA T+17146713818 F+17144948173

sales.usa@apexmedicalcorp.com www.apexmedicalcorp.com Print-2019/05-500 All rights reserved